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效价:  $> 700$  ug/mg

**G-418 Sulfate**

G-418 sulfate is an aminoglycosidic antibiotic, it acts by interfering with protein synthesis and is usually used as an agent for selection of cultured cells expressing a gene that confers resistance to G-418. It is toxic to both prokaryotic and eukaryotic cells. G-418 is white powder, is shipped at room temperature, dry powder should be stored at RT and solution at -20°C.

For research use only. Not for human or drug use.

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分子式:  $C_{20}H_{40}N_4O_{10} \cdot 2H_2SO_4$

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**Kill Curve Assay:**

1. Dissolve G-418 in fully supplemented growth medium without antibiotic, concentration of 5 mg/ml and filter using a 0.22 micron filter.
2. Prepare 6-well cell culture plates by adding G-418 to the growth medium. Concentrations: A range from 100-1,200 µg/ml in 100 µg increments is recommended.
3. Treat cells with trypsin and dilute to a concentration of 4000 cells/ml.
4. Add 100 µl of cell suspension to each well and incubate plates in a humidified CO<sub>2</sub> atmosphere.
5. At 10 to 14 days, aspirate the supernatant and wash the cells with phosphate buffered saline with 0.5% methylene blue and 50% methanol for 20 minutes.
6. Score the plates by calculating percentage of survival by number of individual colonies .
7. Calculate the percentage of survival in the presence of each dilution of G-418 Selective and the percentage of survival in the absence of G-418.
8. Generate a dose response curve by plotting the percentage of survival on the y axis versus G-418 in µg/ml on the x axis for both the sensitive and resistant.

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Sunma Biotech Co., Ltd.