



Sunma

Sofast[®] Transfection Reagent in vivo Protocol

1. Description

Sofast[®] is new generation cationic polymer gene transfection reagent.. Sofast[®] is suitable for transfection of nucleic acids into eukaryotic cells, both in vitro and in vivo gene delivery. Use of Sofast[®] shows high transfection efficiency, low cytotoxicity and simple procedure, it is not cleared by plasma, Sofast[®] shows more advantages in vivo gene delivery compared to viral vector and cationic lipids.

2. Package and storage

Sofast[®] is provided in liquid form at a concentration of 5mg/ml.

Sofast[®] is shipped at room temperature and should be stored at 4°C upon arrival. It is stable for one year at 4°C. Shake gently before use.

3. Transfection procedure (for injection of adult mouse tail vein)

3.1 Dilute 50 µg of DNA into 200 µl of 5% glucose (w/v), mix gently.

3.2 Dilute 7.5 µl Sofast[®] reagent into 200µl of 5% glucose (w/v),mix gently.

3.3 Add 200ul Sofast[®] solution to 200ul DNA solution drop while with vortex.

Notes: The order of mixing two solutions is very important for gene transfection results. Do not reverse the order.

3.4 Incubate for 15-20 minutes at room temperature.

3.5 Inject animals

3.6 Monitor aimed gene expression in desired time. Gene delivery and expression may require 12-48 hours.

Table 1. Recommended amounts of DNA and Sofast[®] according to the route of injection

Animal	Site of injection	Amount of DNA (ug)	Volume of Sofast [®] (ul)	Maximum injection volume (ul)
Adult mouse	Tail vein	50	7-8	200-400
	Portal vein	50	7-8	1000
	Brain ventricle	2.5	0.2-0.3	5
	Heart	50	7-8	200
	Subcutaneous tumor	10	1.2-1.8	100
New born mouse	Brain ventricle	1	0.1-0.15	2

4. Trouble shooting

Problems	suggestions
Too low transfection efficiency	<ol style="list-style-type: none"> 1. Use optimal amount of plasmid 2. Use high quality plasmid, free of RNA (OD_{260/280}>1.8). 3. Optimize the Sofast[®]/DNA ratio: from 1ul Sofast[®]/10ug DNA to 2ul Sofast[®]/10ug DNA.
Animal death	<ol style="list-style-type: none"> 1. Decrease the amount of plasmid (Keep Sofast[®]/DNA ratio constant) 2. Make sure the plasmid is free of endotoxin.